

DemeCAPRONE™

Polyglactone 25 Suture

DemeTECH'S Polyglactone 25 Sutures



DemeTech's fast absorbing Polyglactone 25 is made of a synthetic absorbable monofilament composed of glycolide and ε-caprolactone polymer (PGLA).

DemeTech's Polyglactone 25 is available in either coated or uncoated.

Characteristics of DemeTech's Polyglactone 25 Sutures:

- Very high tensile strength initially
- Excellent handling properties and high pliability
- Very low incidence of infection and trauma
- Smooth tissue passage
- Facile absorption
- High in vivo strength retention



CODES DemeCAPRONE™

No Needle									
Needle Length	Needle Code	Color	Tensile Length	1	2	30	60	60	
			60cm	WGL001	WGL002	WGL003	WGL004	WGL005	
			90cm	WGL006	WGL007	WGL008	WGL009	WGL010	
			120cm	WGL011	WGL012	WGL013	WGL014	WGL015	
Blunt Point 1/2 Circle									
				60					
	WGL016		60cm	WGL016					
				90					
	WGL017		90cm	WGL017					
				120					
	WGL018		120cm	WGL018					
Blunt Point 3/4 Circle									
				60					
	WGL019		60cm	WGL019					
				90					
	WGL020		90cm	WGL020					
				120					
	WGL021		120cm	WGL021					
Precision Curved Standard Point 1/2 Circle									
				60					
	WGL022		60cm	WGL022					
Precision Curved Standard Point 3/4 Circle									
				60			90		
	WGL023		60cm	WGL023			WGL023		
	WGL024		90cm	WGL024			WGL024		
				60			90		
	WGL025		60cm	WGL025			WGL025		
Precision Point Reverse Cutting 1/2 Circle									
				60					
	WGL026		60cm	WGL026					
Precision Point Reverse Cutting 3/4 Circle									
				60			90		
	WGL027		60cm	WGL027			WGL027		
	WGL028		90cm	WGL028			WGL028		
Precision Point Reverse Cutting 3/8 Circle									
				60			90		
	WGL029		60cm	WGL029			WGL029		
	WGL030		90cm	WGL030			WGL030		
				60			90		
	WGL031		60cm	WGL031			WGL031		
	WGL032		90cm	WGL032			WGL032		
				60			90		
	WGL033		60cm	WGL033			WGL033		
	WGL034		90cm	WGL034			WGL034		
Reverse Cutting 1/2 Circle									
				60			90		
	WGL035		60cm	WGL035			WGL035		
	WGL036		90cm	WGL036			WGL036		
				60			90		
	WGL037		60cm	WGL037			WGL037		
	WGL038		90cm	WGL038			WGL038		
				60			90		
	WGL039		60cm	WGL039			WGL039		
	WGL040		90cm	WGL040			WGL040		
Reverse Cutting Straight									
				60			90		
	WGL041		60cm	WGL041			WGL041		
	WGL042		90cm	WGL042			WGL042		
Round Bodied 1/2 Circle									
				60			90		
	WGL043		60cm	WGL043			WGL043		
	WGL044		90cm	WGL044			WGL044		
				60			90		
	WGL045		60cm	WGL045			WGL045		
	WGL046		90cm	WGL046			WGL046		
				60			90		
	WGL047		60cm	WGL047			WGL047		
	WGL048		90cm	WGL048			WGL048		
Round Bodied 1/2 Circle (DOUBLE)									
				60			90		
	WGL049		60cm	WGL049			WGL049		
	WGL050		90cm	WGL050			WGL050		
Round Bodied 1/2 Circle Bevel									
				60			90		
	WGL051		60cm	WGL051			WGL051		
	WGL052		90cm	WGL052			WGL052		
				60			90		
	WGL053		60cm	WGL053			WGL053		
	WGL054		90cm	WGL054			WGL054		
Round Bodied Down Beak 1/2 Circle									
				60			90		
	WGL055		60cm	WGL055			WGL055		
Round Bodied Down Beak 1/2 Circle (DOUBLE)									
				60			90		
	WGL056		60cm	WGL056			WGL056		
	WGL057		90cm	WGL057			WGL057		
Round Bodied 3/4 Circle									
				60					
	WGL058		60cm	WGL058					
				90					
	WGL059		90cm	WGL059					
Round Bodied 3/8 Circle									
				60			90		
	WGL060		60cm	WGL060			WGL060		
	WGL061		90cm	WGL061			WGL061		
Super Cutting 1/2 Circle									
				60			90		
	WGL062		60cm	WGL062			WGL062		
	WGL063		90cm	WGL063			WGL063		