


DemeDAC™

Polyester Suture (Braided without covering)


DemeTECH'S Polyester Sutures
















DemeTech's Polyester (braided without covering) suture is a non-absorbable, sterile, surgical suture composed of Poly (ethylene terephthalate). DemeTech's Polyester (braided without covering) suture elicits a minimal acute inflammatory reaction in tissues, followed by gradual encapsulation of the suture. High quality materials, combined with precise braiding mechanisms, nearly eliminate the occurrence of post-operative suture fragments in the tissue. Polyester (braided without covering) Suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures. Product is intended for long term use over 30 days. This product is not intended to be used on the heart or with the central circulatory and central nervous systems.

Characteristics of DemeTech's Polyester (braided without covering) Sutures:

- As Polyester (braided without covering) is not coated, it allows for extreme ease of use and secure knot tying.
- It provides for minimal breakage, and thus eliminates the need to remove irritating suture fragments postoperatively.
- No significant change in tensile strength known to occur in vivo
- As Polyester (braided without covering) is uncoated, it maintains a higher coefficient of friction as it

Indefinite	Multifilament  Green White	Non Absorbable Suture
Tensile Strength	Material	Absorption Profile

CODES DemeDAC™

Blunt Point 1/2 Circle (DOUBLE)						
Needle Length	Needle Code	Color	Thread Length	5		
65 mm 	DBP-1 ○	White	30 cm	PT850565C16P		
80 mm  <small>*Not to scale</small>	DBP ○	White	30 cm	PT850580C16P		
Conventional Spatula 1/4 Circle (DOUBLE)						
Needle Length	Needle Code	Color	Thread Length	4/0	5/0	6/0
8 mm 	DS-14 ▽	White Green	45 cm	PT1540Q,8E9P	PT1550Q,8E9P	PT1560Q,8E9P PT1460Q,8E9P
Precision Point Reverse Cutting 3/8 Circle						
Needle Length	Needle Code	Color	Thread Length	6/0		
11 mm 	DP-1 ▽	White	45 cm	PT156011F13M		
Needle Length	Needle Code	Color	Thread Length	4/0	5/0	
13 mm 	DP-3 ▽	White	45 cm	PT154013F13M	PT155013F13M	
Reverse Cutting 3/8 Circle						
Needle Length	Needle Code	Color	Thread Length	4/0		
19 mm 	DFS-2 ▽	White	45 cm	PT154019F4P		
Needle Length	Needle Code	Color	Thread Length	3/0		
24 mm 	DFS-1 ▽	White	45 cm	PT153024F4P		
Needle Length	Needle Code	Color	Thread Length	2/0		
26 mm 	DFS ▽	White	45 cm	PT152026F4P		
Reverse Cutting 3/8 Circle (DOUBLE)						
Needle Length	Needle Code	Color	Thread Length	1		
75 mm  <small>*Not to scale</small>	DLR ▽	Green	75 cm	PT240175G4P		
Round Bodied 1/2 Circle						
Needle Length	Needle Code	Color	Thread Length	2/0	3/0	
26 mm 	DSH ●	Green	75 cm	PT242026B0P	PT243026B0P	
Needle Length	Needle Code	Color	Thread Length	0		
36 mm 	DCT-1 ●	Green	75 cm	PT240036B0P		
Round Bodied 1/2 Circle (DOUBLE)						
Needle Length	Needle Code	Color	Thread Length	5		
48 mm 	DCTX ●	White	30 cm 45 cm	PT850548C0P PT150548C0P		
Round Bodied Heavy 1/2 Circle						
Needle Length	Needle Code	Color	Thread Length	5		
36 mm 	DMO-4 ●	White	30 cm	PT850536B1P		