

DemeGUT PLAIN™

Plain Catgut Suture

DemeTECH'S Plain Catgut Sutures



DemeTech's Plain Catgut suture is an absorbable, sterile, surgical suture composed of highly purified connective tissue (mostly collagen) derived from either beef or sheep intestines. Due to an extremely short absorption period, DemeTech's Plain Catgut elicits only minimal reaction in tissues. The percentage of collagen in the catgut suture often determines the quality of the suture. Higher percentages of collagen allow for: superior tensile strength, longer absorption times, and lower reactions in vivo.

Characteristics of DemeTech's Plain Catgut Sutures:

- DemeTech manufactures its Chronic Catgut sutures with the highest quality collagen, maintaining a composition of 97%-98% of pure collagen.
- DemeTech's Plain Catgut is rapidly absorbed.
- Tensile Strength is maintained for about 10 days.
- Absorption is complete within 70 days.
- DemeTech's Plain Catgut is ideal for tissues which heal rapidly and require minimal support.

Tensile strength retained up to 10 days	Monofilament	Complete Absorption 70 Days
Tensile Strength	Material	Absorption Profile

CODES DemeGUT PLAIN™

No Needles								
Needle Length	Needle Code	Color	Thread Length	1	0	20	30	40
		Blue	100 cm					07121000
			120 cm	071801	071800	071820	071830	071840
			150 cm	071901	071900	071920	071930	071940
Conventional Spatula 1/4 Circle (DOUBLE)								
Needle Length	Needle Code	Color	Thread Length	40	50	60		
	07111	Blue	40 cm		07110100P			
			45 cm	07110100P	07110150P	07110200P		
Curved Cutting 1/2 Circle								
Needle Length	Needle Code	Color	Thread Length	30				
	07131	Blue	75 cm		07130100P			
			80 cm	07130100P	07130200P			
Curved Cutting 3/8 Circle (DOUBLE)								
Needle Length	Needle Code	Color	Thread Length					
	07151	Blue	45 cm	07150100P				
			50 cm	07150200P				
Micropoint Curved Spatula 3/8 Circle (DOUBLE)								
Needle Length	Needle Code	Color	Thread Length					
	071602	Blue	45 cm	07160100P				
			50 cm	07160200P				
Micropoint Reverse Cutting 3/8 Circle (DOUBLE)								
Needle Length	Needle Code	Color	Thread Length					
	07171	Blue	45 cm	07170100P				
			50 cm	07170200P				
Precision Cosmetic Standard Prime 3/8 Circle								
Needle Length	Needle Code	Color	Thread Length	50				
	07181	Blue	45 cm	0718010100	0718010100			
			50 cm	0718010100	0718010100			
Precision Point Reverse Cutting 1/2 Circle								
Needle Length	Needle Code	Color	Thread Length					
	07191	Blue	45 cm	0719010100				
			50 cm	0719010100	0719010100			
Precision Point Reverse Cutting 3/8 Circle								
Needle Length	Needle Code	Color	Thread Length					
	07201	Blue	45 cm	0720010100				
			50 cm	0720010100	0720010100			
		Blue	55 cm	0720020100				
			60 cm	0720020100	0720020100			
Reverse Cutting 1/2 Circle								
Needle Length	Needle Code	Color	Thread Length					
	07211	Blue	45 cm	0721010100				
			50 cm	0721010100	0721010100			
		Blue	55 cm	0721020100				
			60 cm	0721020100	0721020100			
Reverse Cutting 3/8 Circle								
Needle Length	Needle Code	Color	Thread Length					
	07221	Blue	45 cm	0722010100				
			50 cm	0722010100	0722010100			
		Blue	55 cm	0722020100				
			60 cm	0722020100	0722020100			
Reverse Cutting Straight								
Needle Length	Needle Code	Color	Thread Length					
	07231	Blue	45 cm	0723010100				
			50 cm	0723010100	0723010100			
		Blue	55 cm	0723020100				
			60 cm	0723020100	0723020100			
Reverse Cutting 3/8 Circle (DOUBLE)								
Needle Length	Needle Code	Color	Thread Length					
	07240	Blue	45 cm	0724010100				
			50 cm	0724020100				
Round Bodied 1/2 Circle								
Needle Length	Needle Code	Color	Thread Length					
	07251	Blue	75 cm	07250100P				
			80 cm	07250200P				
Round Bodied 3/8 Circle								
Needle Length	Needle Code	Color	Thread Length	0	20	30	40	
	07261	Blue	75 cm	07260100P	07260200P	07260300P	07260400P	
			80 cm	07260100P	07260200P	07260300P	07260400P	
Round Bodied 1/2 Circle								
Needle Length	Needle Code	Color	Thread Length	1	0	20	30	
	07271	Blue	75 cm	07270100P	07270200P	07270300P	07270400P	
			80 cm	07270100P	07270200P	07270300P	07270400P	
Round Bodied 3/8 Circle								
Needle Length	Needle Code	Color	Thread Length	1	0	20	30	
	07281	Blue	75 cm	07280100P	07280200P	07280300P	07280400P	
			80 cm	07280100P	07280200P	07280300P	07280400P	
Round Bodied Heavy 1/2 Circle								
Needle Length	Needle Code	Color	Thread Length					
	07290	Blue	75 cm	07290100P				
			80 cm	07290200P				
Round Bodied 3/8 Circle								
Needle Length	Needle Code	Color	Thread Length					
	07301	Blue	75 cm	07300100P				
			80 cm	07300200P				
Taper Cutting 1/2 Circle								
Needle Length	Needle Code	Color	Thread Length					
	07311	Blue	45 cm	07310100P				
			50 cm	07310100P	07310200P	07310300P	07310400P	