


DemeLON MULTI™

Nylon Suture (Multifilament)

DemeTECH'S Nylon Sutures



DemeTech's Nylon suture is a non-absorbable, sterile surgical suture composed of the long chain aliphatic polymers Nylon 6 and Nylon 6.6. DemeTech's Nylon Multifilament sutures are available in either black or colorless, dyed with FDA listed dye(s). This process is usually followed by gradual encapsulation of the suture by fibrous connective tissue. The DemeTECH Nonabsorbable Nylon Suture with and without needles is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures. Product is intended for long term use over 30 days. This product is not intended to be used on the heart or with the central circulatory and central nervous systems.

Characteristics of DemeTech's Nylon Sutures:

- Composed of tightly braided filaments
- Easy passage through tissues due to its low frictional coefficient.
- Minimal tissue reactivity
- Coated to improve handling.
- Excellent elasticity & exceptional for skin closures.
- Small USP sizes are designed specifically for microsurgical and ophthalmic procedures.
- Although not absorbed, progressive hydrolysis may result in gradual loss of tensile strength over time.

Gradual loss of tensile strength over time

Tensile Strength

Multifilament



Black Colorless

Material

Non Absorbable Suture

Absorption Profile

CODES DemeLON MULTI™

No Needle											
Needle Length	Needle Code	Color	Thread Length	4	2	1	0	2/0	3/0	4/0	5/0
		Black	2x20 cm					NM00020			
			2x45 cm					NM00045			
			10x45 cm					NM10045			
			13x45 cm		NM13045	NM13045	NM13045	NM13045	NM13045	NM13045	NM13045
			17x45 cm			NM17045		NM17045			
			2x60 cm	NM20060							
			10x75 cm		NM10075	NM10075	NM10075	NM10075	NM10075	NM10075	
			13x75 cm			NM13075	NM13075	NM13075			
			100 cm	NM00100							
			2x100 cm	NM20100				NM20100			
Curved Cutting 3/8 Circle											
Needle Length	Needle Code	Color	Thread Length	2/0							
26 mm	DSB-25	Black	45 cm	NM162026F2P							
Precision Cosmetic Standard Prime 3/8 Circle											
Needle Length	Needle Code	Color	Thread Length	5/0				6/0			
13 mm	DFC-1	Black	45 cm	NM165013F11M				NM166013F11M			
Reverse Cutting 1/2 Circle											
Needle Length	Needle Code	Color	Thread Length	3/0							
22 mm	DX-1	Black	45 cm	NM163022B4P							
Needle Length	Needle Code	Color	Thread Length	2/0				3/0			
26 mm	DCP-2	Black	45 cm	NM162026B4P				NM163026B4P			
Needle Length	Needle Code	Color	Thread Length	0				2/0			
40 mm	DCP	Black	8x45 cm					NM0x162040B4P			
			90 cm	NM100040B4P							
Reverse Cutting 1/2 Circle											
Needle Length	Needle Code	Color	Thread Length	1							
45 mm	DCF45	Black	90 cm	NM190045B4P							
Reverse Cutting Heavy 1/2 Circle											
Needle Length	Needle Code	Color	Thread Length	1							
22 mm	DOS-4	Black	75 cm	NM260122B5P							
Needle Length	Needle Code	Color	Thread Length	1							
40 mm	DOS-8	Black	3x45 cm	NM0x160140B5P							
Reverse Cutting 3/8 Circle											
Needle Length	Needle Code	Color	Thread Length	2/0				4/0			
26 mm	DF5	Black	45 cm	NM162026F4P				NM164026F4P			
Needle Length	Needle Code	Color	Thread Length	3/0							
30 mm	DF5L	Black	8x45 cm	NM0x162030F4P							
Reverse Cutting Straight											
Needle Length	Needle Code	Color	Thread Length	2/0							
55 mm	DS-55	Black	100 cm	NM162055M4P							
Round Bodied 1/2 Circle											
Needle Length	Needle Code	Color	Thread Length	4/0							
13 mm	DRB-2	Black	45 cm	NM164013B0P							
Needle Length	Needle Code	Color	Thread Length	3/0				4/0			
17 mm	DRB-1	Black	8x45 cm	NM163017B0P				NM164017B0P			
Round Bodied 1/2 Circle											
Needle Length	Needle Code	Color	Thread Length	3/0							
20 mm	DSH-2	Black	75 cm	NM263020B0P							
Needle Length	Needle Code	Color	Thread Length	2/0				3/0			
22 mm	DSH-1	Black	45 cm	NM162022B0P				NM163022B0P			
Needle Length	Needle Code	Color	Thread Length	2/0		3/0		4/0			
26 mm	DSH	Black	45 cm	NM162026B0P		NM163026B0P		NM164026B0P			
			75 cm	NM062026B0P		NM263026B0P					
Needle Length	Needle Code	Color	Thread Length	1		0		2/0			
26 mm	DC3-1	Black	45 cm	NM160126B0P		NM160026B0P		NM162026B0P			
			75 cm	NM060126B0P		NM260026B0P		NM262026B0P			
Needle Length	Needle Code	Color	Thread Length	1				0			
48 mm	DCX	Black	8x45 cm	NM0x160148B0P				NM0x160048B0P			
Round Bodied Heavy 1/2 Circle											
Needle Length	Needle Code	Color	Thread Length	0				2/0			
22 mm	DMO-7	Black	45 cm	NM160022B1P				NM162022B1P			
Needle Length	Needle Code	Color	Thread Length	0				2/0			
26 mm	DMO-6	Black	45 cm	NM160026B1P				NM162026B1P			
Needle Length	Needle Code	Color	Thread Length	1							
30 mm	DMO-5	Black	75 cm	NM260130B1P							
Round Bodied J Needle											
Needle Length	Needle Code	Color	Thread Length	1							
31 mm	DASH-30	Black	75 cm	NM260131J0P							
Taper Cutting Heavy 1/2 Circle											
Needle Length	Needle Code	Color	Thread Length	1							
40 mm	DV-37H	Black	100 cm	NM060140E2P							