


HNM

EMPOWERING INNOVATION



SILICONE GRIPS


Get a sure grip with silicone handles from HNM

Silicone grips from HNM are ergonomically designed to provide a comfortable grip during surgical procedures.

- Excellent resistance to alkaline, saline and acidic solutions.
- Available in different colors.

Custom options:

- connection
- with or without cannulation
- with or without ratcheting
- with or without working end
- with or without quick coupling



MQC = mini quick coupling
 SQC = small quick coupling
 SQQC = square quick coupling
 LQC = large quick coupling
 or custom quick coupling

- handle color
- handle shape

- Non-fading vibrant colors withstand repeated autoclaving.
- Stainless steel core.

Custom designs can be manufactured in a timely and economical manner.



- with or without your company's logo on instrument

SILICONE GRIP

- HNM-3-10-9401 4" black
- HNM-3-10-9402 4" blue
- HNM-3-10-9403 4" red
- HNM-3-10-9404 4" green
- HNM-3-10-9405 4" yellow
- HNM-3-10-9406 4" orange
- HNM-3-10-9407 4" grey

- Connection m6**
- HNM-3-10-9451 4 1/2" black
 - HNM-3-10-9452 4 1/2" blue
 - HNM-3-10-9453 4 1/2" red
 - HNM-3-10-9454 4 1/2" green
 - HNM-3-10-9455 4 1/2" yellow
 - HNM-3-10-9456 4 1/2" orange
 - HNM-3-10-9457 4 1/2" grey



SILICONE GRIP

- Connection m8**
- HNM-3-10-9521 5 1/4" black
 - HNM-3-10-9522 5 1/4" blue
 - HNM-3-10-9523 5 1/4" red
 - HNM-3-10-9524 5 1/4" green
 - HNM-3-10-9525 5 1/4" yellow
 - HNM-3-10-9526 5 1/4" orange
 - HNM-3-10-9527 5 1/4" grey



SILICONE GRIP

- mini quick coupling (MQC)**
- HNM-3-10-9591 5 1/4" black
 - HNM-3-10-9592 5 1/4" blue
 - HNM-3-10-9593 5 1/4" red
 - HNM-3-10-9594 5 1/4" green
 - HNM-3-10-9595 5 1/4" yellow
 - HNM-3-10-9596 5 1/4" orange
 - HNM-3-10-9597 5 1/4" grey



SILICONE GRIP

- small quick coupling (SQC)**
- HNM-3-10-9561 5 3/4" black
 - HNM-3-10-9562 5 3/4" blue
 - HNM-3-10-9563 5 3/4" red
 - HNM-3-10-9564 5 3/4" green
 - HNM-3-10-9565 5 1/4" yellow
 - HNM-3-10-9566 5 3/4" orange
 - HNM-3-10-9567 5 3/4" grey





SILICONE GRIP

- small quick coupling (SOC)*
- HNM-3-10-9571 4 3/4" black
 - HNM-3-10-9572 4 3/4" blue
 - HNM-3-10-9573 4 3/4" red
 - HNM-3-10-9574 4 3/4" green
 - HNM-3-10-9575 4 3/4" yellow
 - HNM-3-10-9576 4 3/4" orange
 - HNM-3-10-9577 4 3/4" grey

cannulated 1.8 mm



SILICONE GRIP

- small quick coupling (SOC)*
- HNM-3-10-9581 5" black
 - HNM-3-10-9582 5" blue
 - HNM-3-10-9583 5" red
 - HNM-3-10-9584 5" green
 - HNM-3-10-9585 5" yellow
 - HNM-3-10-9586 5" orange
 - HNM-3-10-9587 5" grey

with stainless steel turning knob



SILICONE GRIP

- small quick coupling (MQC)*
- HNM-3-10-9531 6" black
 - HNM-3-10-9532 6" blue
 - HNM-3-10-9533 6" red
 - HNM-3-10-9534 6" green
 - HNM-3-10-9535 6" yellow
 - HNM-3-10-9536 6" orange
 - HNM-3-10-9537 6" grey

ratcheting, cannulated 2.3 mm



SILICONE GRIP

- square quick coupling (SQOC)*
- HNM-3-10-9541 6 1/2" black
 - HNM-3-10-9542 6 1/2" blue
 - HNM-3-10-9543 6 1/2" red
 - HNM-3-10-9544 6 1/2" green
 - HNM-3-10-9545 6 1/2" yellow
 - HNM-3-10-9546 6 1/2" orange
 - HNM-3-10-9547 6 1/2" grey

ratcheting, cannulated 3 mm

SILICONE GRIP

large quick coupling (LOC)

- HNM-3-10-9551 6" black
- HNM-3-10-9552 6" blue
- HNM-3-10-9553 6" red
- HNM-3-10-9554 6" green
- HNM-3-10-9555 6" yellow
- HNM-3-10-9556 6" orange
- HNM-3-10-9557 6" grey



ratcheting, cannulated 3.0 mm



SILICONE T-GRIP

connection m6


- HNM-3-10-9611 3 3/4" black
- HNM-3-10-9612 3 3/4" blue
- HNM-3-10-9613 3 3/4" red
- HNM-3-10-9614 3 3/4" green
- HNM-3-10-9615 3 3/4" yellow
- HNM-3-10-9616 3 3/4" orange
- HNM-3-10-9617 3 3/4" grey



SILICONE T-GRIP

connection m8

- HNM-3-10-9621 4 1/2" black
- HNM-3-10-9622 4 1/2" blue
- HNM-3-10-9623 4 1/2" red
- HNM-3-10-9624 4 1/2" green
- HNM-3-10-9625 4 1/2" yellow
- HNM-3-10-9626 4 1/2" orange
- HNM-3-10-9627 4 1/2" grey




SILICONE T-GRIP

connection m8

- HNM-3-10-9641 4 1/2" black
- HNM-3-10-9642 4 1/2" blue
- HNM-3-10-9643 4 1/2" red
- HNM-3-10-9644 4 1/2" green
- HNM-3-10-9645 4 1/2" yellow
- HNM-3-10-9646 4 1/2" orange
- HNM-3-10-9647 4 1/2" grey


cannulated 6.5 mm



SILICONE T-GRIP

- small quick coupling (SQC)*
- HNM-3-10-9701 3" black
 - HNM-3-10-9702 3" blue
 - HNM-3-10-9703 3" red
 - HNM-3-10-9704 3" green
 - HNM-3-10-9705 3" yellow
 - HNM-3-10-9706 3" orange
 - HNM-3-10-9707 3" grey


cannulated 1.8 mm



SILICONE T-GRIP

- square quick coupling (SQQC)*
- HNM-3-10-9751 4 1/2" black
 - HNM-3-10-9752 4 1/2" blue
 - HNM-3-10-9753 4 1/2" red
 - HNM-3-10-9754 4 1/2" green
 - HNM-3-10-9755 4 1/2" yellow
 - HNM-3-10-9756 4 1/2" orange
 - HNM-3-10-9757 4 1/2" grey


ratcheting, cannulated 3.0 mm



SILICONE T-GRIP

- large quick coupling (LQC)*
- HNM-3-10-9721 4 1/2" black
 - HNM-3-10-9722 4 1/2" blue
 - HNM-3-10-9723 4 1/2" red
 - HNM-3-10-9724 4 1/2" green
 - HNM-3-10-9725 4 1/2" yellow
 - HNM-3-10-9726 4 1/2" orange
 - HNM-3-10-9727 4 1/2" grey

ratcheting



SILICONE T-GRIP

- large quick coupling (LQC)*
- HNM-3-10-9741 4 1/2" black
 - HNM-3-10-9742 4 1/2" blue
 - HNM-3-10-9743 4 1/2" red
 - HNM-3-10-9744 4 1/2" green
 - HNM-3-10-9745 4 1/2" yellow
 - HNM-3-10-9746 4 1/2" orange
 - HNM-3-10-9747 4 1/2" grey

ratcheting, cannulated 3.0 mm

SILICONE TEARDROP GRIP


- connection m8
- HNM-3-10-9821 3 1/2" black
 - HNM-3-10-9822 3 1/2" blue
 - HNM-3-10-9823 3 1/2" red
 - HNM-3-10-9824 3 1/2" green
 - HNM-3-10-9825 3 1/2" yellow
 - HNM-3-10-9826 3 1/2" orange
 - HNM-3-10-9827 3 1/2" grey



SILICONE TEARDROP GRIP

- square quick coupling (SQC)
- HNM-3-10-9841 5" black
 - HNM-3-10-9842 5" blue
 - HNM-3-10-9843 5" red
 - HNM-3-10-9844 5" green
 - HNM-3-10-9845 5" yellow
 - HNM-3-10-9846 5" orange
 - HNM-3-10-9847 5" grey

cannulated 6.5 mm



SILICONE WRENCH

- HNM-3-10-9901 10 1/2" black
- HNM-3-10-9902 10 1/2" blue
- HNM-3-10-9903 10 1/2" red
- HNM-3-10-9904 10 1/2" green
- HNM-3-10-9905 10 1/2" yellow
- HNM-3-10-9906 10 1/2" orange
- HNM-3-10-9907 10 1/2" grey

ratcheting ring, SW 17 mm



HNM MEDICAL
TOLL FREE 1.866.291.8498
FAX 1.888.608.9788
WWW.HNMMEDICAL.COM
E-MAIL: INFO@HNMMEDICAL.COM